

Hydro-Bandaïd



Nakyung Lee, Qingyuan Liu, Ayako Ohoka, Arushi Sardana, Aaron
Schwartz-Duval, Veronika Sowers

Melanoma



What it is

Most dangerous form of skin cancer originating from the pigment-producing melanocytes in the bottom layer of the epidermis

Causes


Unrepaired DNA damage to skin cells (usually from **ultraviolet radiation**) → genetic mutations → skin cells multiply rapidly

Diagnosis

ABCDEs of melanoma - asymmetry, border, color, diameter, evolving
Biopsy

Stages and Current Treatments

Stage 0 (in situ)	Melanoma only in top surface layer of epidermis Has not metastasized	Wide-local excision surgery
Stage 1	Melanoma less than 2 mm thick Has not metastasized	Wide-local excision surgery
Stage 2	Melanoma between 1 mm and around ~4 mm thick Has not metastasized	Wide-local excision surgery



- Skin cancer and a small margin of healthy tissue around it is cut out
- Edges of wound are sutured
- Tissue sent to pathologist for processing and margin evaluation
- Recovery rate depends on size and site of cancer

Significance

Melanoma

- ~73,870 diagnosed in 2015
 - ~9,940 deaths in 2015
- Incidence of melanoma has been increasing for the past 30 years
- Recurrence occurs 10 or more years after initial treatment in more than 1 in 20 patients

Surgical Site Infection

- 3 in every 100 patients who have surgery develop a surgical site infection
 - Most common health-associated infection
 - ~157,500 occurrences a year
- 98,987 deaths by healthcare associated infection in 2002

Our Approach: Hydro-Bandaid


- Focusing on patient recovery post wide-local excision surgery
- Sticky degradable hydrogel patch with microneedle for diffusion into skin
- Layered assembly of alternating hydrogels with environmentally sensitive structure
- Offers controlled release of therapeutics
 - Imiquimod**- prevents recurrence
 - Antibiotics**- prevents infection
 - Interferon**- boosts immune system
 - Anticoagulants**- prevents blood-clotting



Temperature and pH Dependence


poly(n-isopropylacrylamide) – Temperature Sensitive

- Modifiable with cellulose
- Skin~34°C, Internal Body Temp=37°C



poly(acrylic acid) – pH Sensitive

- Modifiable with cellulose
- Skin~5.5, Blood=7.4, Tumor=6.0-6.8




Layered Assembly versus Hybrid Gel

Hybrid Gel

- Uneven dissolution/ release of therapeutic
- Non-uniform diffusion gradient of heat and protons



Layered Assembly

- Each layer provides “insulation” against the affector of the next layer
- Superior control for directional release through two mechanisms
- Separate tunable distinct modes of release for each layer





Thick versus Thin Layered Assembly



Thin layers will result in a more stable release of therapeutic, but the “insulation effect” will be diminished



Thick layers will result in a more pulsatile release of therapeutic ☐ controlled therapeutic release




Fabrication / Methods



- Hydrogel can be prepared in large-scale sheets ☐ Mass production possible
- Cut into desired size using x-y directional laser cutter

Fabrication / Methods



- Layers of hydrogel and drug stacked up
- Hydrogel with drug incorporated in goes on top of the metal micro-needle

- Painless microneedles
- Stainless steel

Testing: Material

Physical Characterization

- Strip extensimetry
- Compression
 - Tensile modulus
- SEM

Drug Release


- Beaker containing 37°C pH 5.5 citrate buffer covered by dialysis membrane
- Read drug concentration in the solution overtime with UV/V is over various timepoints



Testing: In vitro

Direct Contact Test:


- Place hydrogels in well plate containing PCS-201-012 human dermal epithelial cells



Live/Dead Assay

Elution:

- Place hydrogels in culture media in a well plate, then remove media after 24 hours
- Take media to culture PCS-201-012 cells



Cell assays:
MTT, LDH, WST





Testing: In vivo

Inject C32 melanoma cells into hairless guinea pig (HGP), remove most of tumor




Adhesion of Hydro-Bandaïd to excision site

- Drug release study from examining concentration in blood
- Histopathology
 - Inflammation
 - Vascularization
 - Fibrosis



skin
sample



Summary

- Wide-local excision surgery most common treatment for melanoma
- Recovery can be extensive and long
- Hydro-Bandaïd offers controlled release of therapeutics to accelerate recovery and prevent recurrence at primary site
- Future applications: controlled release of antibiotics and/or sterilization agents
 - Topical postoperative infection control
 - Steroid treatment
 - Nicotine withdrawal remedy